

# Eliminate Downtime and Calibrate In-House

Valid until 31 March 2025

Regular calibration and verification of your test and measurement equipment helps keep your world up, running and safe. Calibrate in-house with ease to help save you money and downtime with these latest promotions.

## Temperature Calibration

### Fluke Calibration 9100S Handheld Dry-Well



World's smallest, lightest and most portable dry-well.

- Range from 35°C to 375°C
- Accuracy to  $\pm 0.25^\circ\text{C}$ , stability of  $\pm 0.07^\circ\text{C}$  at 50°C
- RS-232 interface

### Fluke Calibration 9102s/9009 Field Metrology Wells



#### 9102S

- Ranges from  $-10^\circ\text{C}$  to  $122^\circ\text{C}$  in one unit

#### 9009

- Ranges from  $-15^\circ\text{C}$  to  $350^\circ\text{C}$  in one unit
- Two wells in each block for simultaneous comparisons

#### Handheld Dry-Well Offer



Receive FREE YETI Roadie 24L Hard Cooler and YETI Rambler 769mL Drink Bottle valued at \$410.00

#### Dual-Well Offer



Receive FREE Fluke 1551A-9 Stik Thermometer & Temperature Calibrator valued at \$1,704.00

### Fluke Calibration 9142/9143/9144 Field Metrology Wells



Lightweight and portable dry wells. Three models available to cover the wide range of  $-25^\circ\text{C}$  to  $660^\circ\text{C}$ .

#### Purchase the Base Unit

Process Electronics version provides built-in thermometer readout that measures PRTs, RTDs, thermocouples, and 4-20mA transmitters.

### Fluke Calibration 6109A/7109A Portable Calibration Baths



Portable calibration baths with four times more calibration throughput and twice the accuracy of Micro-Baths and dry-block calibrators.

- Calibrate up to four tri-clamp sanitary sensors at the same time

#### Field Metrology Wells Offer

Buy a Base Unit and upgrade to Process Electronics for FREE UPGRADE valued up to \$3,323.00

#### Portable Bath Offer

Buy a Base Unit and upgrade to Process Electronics for FREE UPGRADE valued up to \$3,284.00

# Eliminate Downtime and Calibrate In-House

## Electrical Calibration

### Fluke Calibration 1586A Super-DAQ Precision Temperature Scanner and DAQ-STAQ Multiplexer



Best-in-class temperature measurement accuracy and up to 40 isolated input channels.

- PRTs:  $\pm 0.005^{\circ}\text{C}$ ; Thermocouples:  $\pm 0.29^{\circ}\text{C}$ ; Thermistors:  $\pm 0.002^{\circ}\text{C}$
- Scan speed of up to 10 channels per second

### Fluke 8808A Digital Multimeters



Precision and versatility to eliminate production mistakes and handle demanding measurements.

- 5.5 digit resolution
- Basic V DC accuracy of 0.015%
- Dual display

#### Data Acquisition Calibration Scanner Offer



Receive **FREE** Fluke 1586-CASE Carrying Case **or** Fluke 1586/DS-CASE Carrying Case **valued up to \$2,295.00**

#### Precision Multimeters Offer



Receive **FREE** Fluke Electronic Test Probes **valued at \$195.00**

## Pressure Calibration


### Fluke Calibration 700HPPK Pneumatic Test Pump Kit



Rugged and portable way to generate pressure in the field quickly, safely and easily.

- Generates and adjusts pneumatic pressure up to 3000 psi (21 MPa)

### Fluke Calibration P5514B Hydraulic Pressure Comparator



Pump for calibrating gauges, switches and transducers.

- Generates pressures to 10 000 psi (70 MPa)

#### Pneumatic Test Pump Kit Offer



Receive **FREE** Fluke 700G29 Pressure Gauge Calibrator **valued at \$1,730.00**

#### Hydraulic Pressure Comparator Offer



Receive **FREE** Fluke 2700G Pressure Gauge Calibrator **valued at \$2,094.00**

For the full list of models applicable to this promotion, please reach out to your distributor.

# Efficient, accurate calibration

## One solution for every need

### Electrical Calibration and Software

#### Fluke Calibration 5560A, 5550A and 5540A Multi-Product Calibrators



Make calibration simpler and do more with a single calibration tool.

- High accuracy across multiple measurement ranges
- Visual Connection Management™ guides calibration setups with lighted terminals
- 30 A continuous current output
- Now available: oscilloscope calibration options up to 2 GHz

Features	5560A	5550A	5540A
4:1 test uncertainty ratio for most popular 6.5-digit multimeters	•	-	-
4:1 test uncertainty ratio for most popular 5.5-digit multimeters	•	•	-
4:1 test uncertainty ratio for most popular 4.5-digit multimeters	•	•	•

Eligible models include: 5560A, 5550A, 5540A, 5560A/600M, 5560A/1G, 5560A/2G, 5550A/600M, 5550A/1G and 5540A/600M.

### Software Offers

#### Option 1



SAVE  
**50% OFF**

**MET/CAL™**  
**MET/TEAM™**

#### MET/CAL and MET/TEAM Offer

Buy a Fluke Calibration 55xxA Calibrator and **SAVE 50% OFF** MET/CAL and MET/TEAM Licenses

#### Option 2



SAVE  
**35% OFF**

**MET/SUPPORT™ Gold**

#### MET/SUPPORT GOLD Offer

Buy a Fluke Calibration 55xxA Calibrator and **SAVE 35% OFF** your purchase of MET/SUPPORT GOLD Software New or Renewal Work Stations

#### How to redeem your FREE gift:

1. Purchase a Fluke test tool within this promotion between **4 November - 31 March 2025**.
2. Visit [fluke.co/fcal24](http://fluke.co/fcal24) to redeem your bonus



Scan to redeem:  
[fluke.co/fcal24](http://fluke.co/fcal24)

**\*Terms and Conditions:** Sponsored by Fluke Australia Pty Ltd. Gift with purchase is available for all purchases of a "Qualifying Fluke Purchase" made from an authorised participating Fluke distributor in Australia and New Zealand between **4 November and 31 March 2025** or until stocks last. Gift offerings are "up to" the listed dollar value. Prices are in \$AUD and are Fluke Australia's Recommended Trade Price inclusive of GST. Discounts are applied off the list price. Offers cannot be combined with any other discounts, promotions or other specials offers from Fluke Australia or the participating authorised distributor. Offers are available to customers only. Distributors and/or branches cannot redeem the free items. Fluke reserves the right to substitute a gift of equivalent value or reject the submission if it doesn't meet the redemption requirements. All trademarks are the property of their respective owners. Specifications subject to change without notice. Allow 6-12 weeks for gift delivery. Additional terms and conditions apply.