

Eliminate Downtime and Calibrate In-House

Valid until 31 December 2023

Regular calibration and verification of your test and measurement equipment helps keep your world up, running and safe. Calibrate in-house with ease to help save you money and downtime with these latest promotions.

Temperature Calibration

Fluke Calibration 9100S Handheld Dry-Well



World's smallest, lightest and most portable dry-well

- Ranges from 35°C to 375°C
- Accuracy to ± 0.25 °C, stability of ± 0.07 °C at 50°C
- RS-232 interface

Handheld Dry-Well Offer



Receive a FREE
Fluke Calibration 9300
Carrying Case
valued at \$433.00

Fluke Calibration 9102s/9009 Field Metrology Wells



9102S

- Ranges from -10°C to 122°C in one unit
- 9009**
- Ranges from -15°C to 350°C in one unit
 - Two wells in each block for simultaneous comparisons

Dual-Well Offer



Receive a choice of two FREE
Fluke Calibration Inserts
valued at \$749.00

Fluke Calibration 9142/9143/9144 Field Metrology Wells



Lightweight and portable dry wells

- Three models available to cover the wide range of -25°C to 660°C
- The optional process version provides a built-in two-channel thermometer readout that measures PRTs, RTDs, thermocouples, and 4-20mA transmitters.

Field Metrology Wells Offer



Receive a FREE 1551A-12 or
1552-12 Stik Thermometer &
Temperature Calibrator
Value up to \$2,366.00

Fluke Calibration 6109A/7109A Portable Calibration Baths



Portable calibration baths with four times more calibration throughput and twice the accuracy of Micro-Baths and dry-block calibrators

- Calibrate up to four tri-clamp sanitary sensors at the same time

Portable Bath Offer



Receive a FREE 1551A-12 or
1552-12 Stik Thermometer &
Temperature Calibrator
Value up to \$2,366.00

How to redeem your FREE gift:

1. Purchase an eligible Fluke test tool between 11 SEP - 31 DEC 2023.
2. Visit [fluke.co/calih23](https://www.fluke.com/calih23) to redeem.



Scan QR code for more info
➔ [fluke.co/calih23](https://www.fluke.com/calih23)

Terms and Conditions: Sponsored by Fluke Australia Pty Ltd. Gift with purchase is available for all purchases of a "Qualifying Fluke Purchase" made from an authorised participating Fluke distributor in Australia and New Zealand between 11 September - 31 December 2023 or until stocks last. Redemption gift will be provided by your Fluke distributor. Gift offerings are "up to" the listed dollar value. Prices are in \$AUD and are Fluke Australia's Recommended Trade Price inclusive of GST. Customers are limited to two redemption claims only. Offers cannot be combined with any other discounts, promotions or other specials offers from Fluke Australia or the participating authorised distributor. Fluke reserves the right to substitute a gift of equivalent value or reject the submission if it doesn't meet the redemption requirements. All trademarks are the property of their respective owners. Specifications subject to change without notice. Allow 12 weeks for gift delivery. Additional terms and conditions apply.