

Fiber Kit #1

High performance solution which accelerates every step of the fiber certification process, comprising a Quad Optical Loss Test Set.

Save \$6,416

Item no. 4955805 |
Part no. CFP2-100-Q/GLD INT

+Gold Support



Fiber Kit #2

High performance solution comprising a Quad Optical Time Domain Reflectometer (OTDR) and fiber inspection scope for troubleshooting and certification.

Save \$5,219

Item no. 4955462 |
Part no. OFP2-100QI/GLD INT

+Gold Support



FIBER



Fiber Kit #3

Enhance your existing DSX Copper Series (DSX2-5000/DSX2-8000) and deliver industry-leading efficiency for optical loss certification across both multimode and single mode fibers.

Save \$4,425

Item no. 4334851 |
Part no. CFP-Q-Add



COPPER



FIBER



GOLD SUPPORT

Premium Care – Gold Support

Premium Care is our comprehensive maintenance program for Fluke Network's products.

It provides unmatched support for customers wanting to minimize business downtime and ensure a high return on investment. The service agreement includes annual calibration (where relevant), replacement accessories, repairs with priority service, and access to loaner units during repairs. Gold support is included for select discounted products or can be purchased separately. Contact your Fluke Networks representative today to find out more.

Read more:

www.flukenetworks.com/content/gold-product-support

*Terms and conditions apply. The discount offered under this promotion is valid on purchases of qualifying products made from Fluke's authorized distributors only. Our authorized distributors are free to provide additional price discounts and, thus, higher end-customer savings. Promotion sponsor: Fluke Europe BV. For more information, call: +44-(0)20-79420721, +31-(0)20-7132984, or email: info@flukenetworks.com

